ACADEMY-2 CAMPUS-2 INSTRUCTIONS

www.konus.com

EN

**BIOLOGICAL MONOCULAR MICROSCOPE
ACADEMY-2
BIOLOGICAL BINOCULAR MICROSCOPE
CAMPUS-2**
TECHNICAL CHARACTERISTICS

- Biological monocular microscope with pivoting head a 360° (Binocular head for mod. Campus)
- Double focusing device (macrometric and micrometric)
- N.A. 1.2 condenser with iris diaphragm and filter holder
- 10X eyepiece
- Illumination with Led lamp 0.5W and external transformer
- Four achromatic objectives (4x; 10x; 40x retractable and 100x immersion)
- Magnification in equipment: 40x; 100x; 400x and 1000 times
- Movable stage
- Coloured filters
- Anti Dust cover
- A/C adapter
- Rechargeable batteries and battery charger included
- Triple illumination system: works with A/C adapter, rechargeable batteries or standard AAA batteries.

COMPONENTS (SEE IMAGE)

1. Eyepiece
2. Head
3. Stand
4. Base
5. Rotating turret objective holder
6. Objectives
7. Movable stage
8. Specimen table
9. Condenser with iris diaphragm and filter holder
10. Halogen lamp (internal)
11. On/off switch
12. Macrometric focusing
13. Micrometric focusing (fine regulation)

MOUNTING

Open the packaging and take out all the components.

Position the base (4) with the stand (3) and head (2) on a dry and stable table.

- Take the head (2) and insert it on the upper part of the stand (3)
- Take the 10x eyepiece (1) and after taking off the caps, insert it into the head (2)
- Take out the objectives (6) from their casing and screw them onto the turret (5), turn the turret until the 4x objective (the shortest) is directly above the specimen stage (8).
- Take the transformer and plug it into an electric socket, (make sure that the data on the label of the transformer matches that of your electrical supply)
- Insert the jack plug of the transformer into the socket of the base (4) of the instrument.
- Turn on the light with the two way switch (11). Now your

microscope is ready for use.

CHANGING THE BULB

This microscope is supplied with a Led bulb, in case it needs replacing, take the following steps:

1. Remove the jack plug from the microscope.
2. Turn the micrometric focus knob (13) lifting the specimen table (8) until the limit.
3. Move the condenser (9) as high as possible
4. Unscrew the halogen lamp holder
5. Replace the lamp with a new one
6. Screw on back the lamp holder.

HIGH ENLARGEMENTS:

The objectives from 40 to 100x are much more longer than the objectives of 4x and 10x, and so, rotating the objective turret can touch the specimen slide (8), in this case, it is necessary to first lower the specimen stage, so as not to create any damage.

You will notice that during focusing of the 40x and 100x objectives, they will first touch the cover slide and then the terminal section will enter the inside of the objective. This is absolutely necessary for high enlargements and they have been prepared for such operations with a special spring.

Be very careful not to put too much pressure onto the focusing knob as the specimen slide might be cracked by the force.

The 100x objective supplies a high number of enlargements, and in these conditions, the air present around the objective can highly worsen the quality of the image. To take care of this problem, an immersion oil may be used; using one or two drops of the oil on top of the cover slide, so as not to have any air between the objective and the specimen during observation.

During the observation of very transparent samples, the light might be too intense to be able to see all the details. To resolve this problem, the condenser must be adjusted (9) with an iris diaphragm or coloured filters.

Highering up the condenser and/or opening the diaphragm means letting in more light, while lowering the eyepiece and/or closing the diaphragm means decreasing the light. For particular specimens it is possible to use appropriate dyes to bring out the details.

MAINTENANCE:

The instrument does not present itself with particular maintenance problems, it is important to follow this simple rules:

- Do not temper with the instrument, especially the electrical wiring and the optical parts.
- Protect the instrument from dust, acids, shocks and/or impacts, water and/or humidity. Do not place the transformer within reach of children.

The eventual vision by children of the microscope should be taken under the supervision of adults.

OPTIONAL ACCESSORIES:

For these microscopes are available several spare accessories as: eyepieces, photo adapter, objectives and prepared slides, please check on our website: www.konus.com

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www.konus.com - ©KONUS 2017 all right reserved