Levenhuk 450 Series Biological Microscopes


450M monocular | 450B binocular | 450T trinocular | D450L digital

- ⓔⓃ User Manual
- ⓑⓖ Ръководство за потребителя
- ⒸⓏ Návod k použití
- ⓓⓔ Bedienungsanleitung
- ⓔⓢ Guía del usuario
- ⓗⓤ Használati útmutató
- ⓔⓣ Guida all' utilizzo
- ⓖⓗ Instrukcja obsługi
- ⓖⓣ Manual do usuário
- ⓗⓤ Инструкция по эксплуатации
- ⓣⓗ Kullanım kılavuzu

*Наслади се отблизо
Radost zaostřít
Zoom ran und hab Fun!
Amplie y disfrute
Kellemes nagyítást!
Ingrandisci il divertimento
Dê um zoom na sua emoção
Radość przybliżania
Dê um zoom na sua emoção
Приближает с удовольствием
Yakınlaştırın ve Keyfini Çıkarın*

levenhuk
Zoom&Joy

Levenhuk 450M

Levenhuk 450B

Levenhuk 450T

Levenhuk D450L

1. Levenhuk 450 Series Microscopes

EN

1. Eyepiece
2. Eyepiece head
3. Arm
4. Revolving nosepiece
5. Objective
6. Specimen holder
7. Coarse focusing knob
8. Fine focusing knob
9. Stage moving knob (back/forth)
10. Stage moving knob (right/left)
11. Brightness adjustment ring
12. Power button
13. Collector
14. Aperture diaphragm adjustment rod
15. Condenser
16. Stage
17. Locking screw
18. Diopter adjustment ring
19. Trinocular tube
20. LCD screen

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Caution! Please remember that mains voltage in most European countries is 220–240V. If you want to use your device in a country with a different mains voltage standard, remember that use of a converter is absolutely necessary.

General Information

Read the user manual carefully before you start working with a microscope. Levenhuk 450 Series biological microscopes are safe for health, life and property of the consumer and the environment when properly used, and meets the requirements of international standards. These microscopes are designed for observing transparent objects in the transmitted light using the bright field method. These microscopes are shipped with high-quality achromatic objective lenses with 4x, 10x, 40x, and 100x magnifications that provide for sharp and bright images. To obtain high-quality images in microscopy, microscopes use an iris diaphragm. It allows for the control of brightness and contrast of the image, as well as serves to improve depth of field.

All models may be found in schools and research institutes, at home and in manufacturing plants. They perform well in clinical research and tests; teaching demonstrations; bacterioscopy and cytology in medical and health establishments, laboratories, universities, etc.

Microscope assembly

- Unpack the microscope carefully and place it on a flat surface.
- Inspect all the items included in the kit and define their purpose.
- Remove the plastic bags and dustproof cover of the eyepiece head.
- Install the monocular, binocular, trinocular head or LCD screen and tighten the locking screw.
- Insert the eyepieces into the eyepiece tubes. Rotate the eyepieces, making sure they are tightly seated in the tubes. Set the diopter adjustment ring to zero.
- Lower the stage. Insert the objectives into the sockets of the revolving nosepiece in increasing order of magnification.
- Raise the stage. Insert the Abbe condenser into the condenser holder. Secure the condenser in the holder with the screw.
- Open the battery compartment cover (located underneath), position the batteries according to the polarity marking on the battery compartment, reinsert the battery compartment cover and turn on the microscope. **These microscopes use rechargeable batteries only!** Before using the microscope, it is necessary to charge the batteries. Connect the power cable to the microscope and the DC adapter via a USB Type-C plug and connect it to the AC power supply to operate from AC mains and charge the device.

Use

Please refer to the fig. 1 (monocular model, binocular model, trinocular model, and model with LCD screen).

- Turn on the power and the illumination will turn on. Set the brightness at approximately 70%. Place the specimen on the stage. Make sure the clips hold the specimen firmly in place.
- Start your observations with the lowest magnification objective (4x) and select a specimen segment for detailed research. Then move the specimen to center the selected segment in the field of view, to make sure it keeps centered when the objective is changed to a more powerful one. Once the segment is selected, you should center its image in the microscope's field of view as precisely as possible. Otherwise, the desired segment might fail to center in the field of view of the higher power objective.
- Now you can switch to a more powerful objective by rotating the revolving nosepiece. Adjust the image focus, the iris diaphragm and the illumination brightness if required.
- The diaphragm allows setting a specific aperture to work with brightfield objectives. Iris diaphragm should be adjusted in accordance with numerical aperture of a lens. Close the aperture diaphragm to 1/3 of the objective exit pupil using the adjustment knob. This value will be slightly less than the objective aperture. If the adjustment knob is turned to the right, the iris diaphragm is fully open. **Note:** the diaphragm is not intended for adjusting the brightness of the illumination. To adjust the brightness, use the brightness adjustment ring.
- To change the objective magnification, turn the revolving nosepiece until it clicks.
- Make sure that the lens does not touch the sample when adjusting the focus: rotate the coarse focus knob until the specimen is about 3.175mm (1/8") away from the lens.

- Looking into the eyepiece, slowly turn the coarse focusing knob until you can see the image of the observed specimen. Turn the fine focusing knob to make the image clear. The fine focusing mechanism allows you to focus on the observed specimen while using high magnifications. When observing a large number of the same type specimens, the coarse focus locking mechanism can be used. This allows for easy specimen changes without losing the set focus. This is especially convenient when the specimens have similar thickness and optical properties.
- The condenser is height adjustable. Move the condenser upward when the magnification of objective lens increases.
- When you use an oil immersion objective, the space between the specimen and the lens should be filled with immersion oil. Raise the condenser to the upper position and put a drop of oil on the lens and on the cover glass of the specimen. Move the revolving nosepiece back and forth to get rid of air bubbles in the oil. Then, fix the objective lens in its working position. Make sure that oil fills the entire space between the objective lens and the specimen. After use, wipe the lens dry.
- To place the specimen in the field of view, move the stage horizontally back and forth or left and right using the stage moving knobs.
- Look through the right eyepiece with your right eye and adjust the image sharpness with the coarse and fine focusing knobs. Binocular and trinocular models: Then, look through the left eyepiece with your left eye and rotate the diopter adjustment ring to equalize the difference between your left and right eyes vision. Remember the setting for future use. Adjust the distance between the eyepiece tubes of the binocular head so that the image merges into a single circle.
- Trinocular model: you can install the digital camera to the third eyepiece tube.

LCD screen (Levenhuk D450L)

Getting started

- Connect the power cable to the power input on the LCD screen and connect it to the network. The red LED will light up. Turn the LCD screen on with the power on/off button; the LED light will change to green.
- Insert the SD card (included) in the SD slot to view the image and save the resulting photos and videos.

Menu and functions

- Use the buttons on the back of the LCD screen to control the camera functions and parameters. After modifying the parameters, exit the interface to save data.
- Press the Menu button to call up the menu. Use $\uparrow\downarrow$ to select the functions. Use \rightarrow to enter the submenu. Press the Menu button once again to save the modified parameters and exit the interface.

Photo and video


- **Photo.** Press the Snap button on the LCD panel to take a photo. The image will be saved on the SD card.
- **Video.** Before recording a video, make sure that the SD card is formatted with the FAT32 file system and has enough free space. You cannot take photos while recording a video. Enter the menu and select "Record Video" to record a video.

Basic settings

- **White Balance.** In this parameter, you can adjust the color temperature. The default option is "Auto White Balance".
- **Exposure.** In "Auto Exposure", you can set the exposure value. In "Manual Exposure", you can adjust the "time" parameter.
- **Color Adjustment.** In this parameter, you can adjust saturation, hue, brightness, and contrast.
- **Monochrome.** By selecting this parameter, you will see the image in different tones of a single color (e.g. shades of gray).
- **Flip.** This feature allows you to flip the image horizontally or vertically.

Measurements

1. **Cross Line.** You can select a cross line in four colors, turn the display on or off, and adjust the position of the center point of the cross line.
2. **Calibration.** The camera is calibrated by default, but you may need to recalibrate the camera.
 - Put the calibration slide on the stage so that the scale is placed horizontally on the screen. Adjust the focus as clearly as possible.
 - Enter the Calibration menu. Drag a segment of the calibration ruler so that its ends are next to the marks on the physical scale of the calibration slide. For maximum accuracy, try to capture as many scales as possible.

- Example: The image above uses a ruler with 0.01 mm (10µm) divisions under a 10x lens. The camera parameters are set to: "magnification 10x", "units µm", and "length 40".
 - Each objective lens must be calibrated to take measurements at different magnifications.
 - After adjusting the parameters, exit the calibration menu to complete the procedure.
3. **Recalibration.** The calibration process has to be repeated if the magnification or focus has been changed during observations. Enter the Measurement menu. Select the magnification, adjust the start and end point, and then the value of the measured segment will be displayed.

Specifications

	450M	450B	450T	D450L
Magnification, x	40–1000			
Tube length, mm	160			
Microscope head	monocular 360° rotatable; 30° inclined; eyepiece diameter: 23.2mm	binocular Siedentopf type, 360° rotatable; 30° inclined; interpupillary distance: 48–75mm eyepiece diameter: 23.2mm	trinocular Siedentopf type, 360° rotatable; 30° inclined; interpupillary distance: 48–75mm eyepiece diameter: 23.2mm	7" color LCD screen, IPS, 2 Mpx resolution: 1824x1068px
Eyepieces, x/field, mm	10x/18, diopter adjustment: ±5D			–
Revolving nosepiece	4 objectives			
Optical design	achromatic objectives, parfocal distance: 45mm			
Objectives, x/aperture	4x/0.10; 10x/0.25; 40x/0.65; 100xs/1.25 (oil)			
Mechanical stage	double layer mechanical stage, stage size: 110x120mm, moving range: 70x30mm			
Condenser	Abbe condenser, N.A. 1.25, with adjustable iris diaphragm			
Focusing mechanism	coaxial coarse & fine focusing knobs			
Light source	LED 1W			
Power supply	100–240V via 5V/1A adapter, or 3 AA rechargeable batteries			
Operating temperature range, °C	0... +40			
Operating humidity range, %	20... 85			

LCD screen specifications (Levenhuk D450L)

Megapixels	2
Max. resolution (still images)	1824x1068px
Frame rate	30fps
Sensor	1/2.8"
Mounting location	microscope's arm
Pixel size	2.9x2.9µm
Sensitivity	2000mV
Image format	*.jpg
Video format	*.mp4
Exposure time	0.1–33ms
White balance	auto/once/manual
Exposure control	auto/manual
Storage	4GB SD card
Output	SD slot, USB2.0 output connection to PC (optional)
Power supply	5V, 1A via AC adapter
Software	Embedded
Operating system	Windows 11 and below

The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

Care and maintenance

Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS. Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions. After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection. Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center. Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws. Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics. Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush. Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity. Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains. If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope. Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures. When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may heat up during use, and might become a fire hazard. Always unplug the microscope from a power source before opening the base or changing the illumination lamp. Regardless of the lamp type (halogen or incandescent), give it some time to cool down before trying to change it, and always change it to a lamp of the same type. Always use the power supply with the proper voltage, i.e. indicated in the specifications of your new microscope. Plugging the instrument into a different power outlet may damage the electric circuitry of the microscope, burn out the lamp, or even cause a short circuit. Children should use the device under adult supervision only. **Seek medical advice immediately if a small part or a battery is swallowed.**

Battery safety instructions

The device is equipped with a rechargeable lithium battery. This avoids frequent battery replacement. Always switch the device off when not in use. If the battery charge is low, please recharge the device in time to ensure that the measurements are accurate. Do not overheat the battery. Do not discharge the battery completely.

Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars, and other optical products, except for their accessories, carry a **lifetime warranty** against defects in materials and workmanship. A lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from the date of purchase. The warranty entitles you to free repair or replacement of your Levenhuk product in any country where a Levenhuk representative office is located, provided that all terms and conditions of the warranty are met.

For further details, please visit our web site: www.levenhuk.com/warranty

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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Levenhuk prepared slides sets

The Levenhuk slides sets include thoroughly prepared specimens of various biomaterials as well as blank slides and cover glasses for you to make your own specimens. Prepared slides, blank slides and cover slips can be used with any microscope model.

Levenhuk Inc. (USA): 928 E 124th Ave. Ste D, Tampa, FL 33612,
USA, +1 813 468-3001, contact_us@levenhuk.com
Levenhuk Optics s.r.o. (Europe): V Chotejně 700/7, 102 00 Prague 102,
Czech Republic, +420 737-004-919, sales-info@levenhuk.cz
Levenhuk® is registered trademark of Levenhuk, Inc.
© 2006–2024 Levenhuk, Inc. All rights reserved.
www.levenhuk.com
20241112

levenhuk
Zoom&Joy